THUASNE

LYMPHATREX®

de Gebrauchsanweisung - Kompressionsbekleidung zur Ödemtherapie
en Instructions for use - Compression garments for oedema therapy
fr Consignes d'utilisation - Vêtement de compression destiné au traitement des oedèmes
it Istruzioni per l'uso - Indumenti di compressione per la terapia degli edemi
es Instrucciones de utilización - Prenda de compresión para la terapia de edemas
pt Instruções de uso - Vestuário de compressão para terapia de edemas
nl Gebruiksaanwijzing - Compressiekleding voor de oedeemtherapie
da Brugsanvisning - Kompressionstøj til ødembehandling
sv Bruksanvisning - Kompressionsplagg för ödembehandling
fi Käyttöohje - Kompressiovaatteet turvotuksen hoitoon
no Bruksanvisning - Kompresjonsbekledning til ødemterapi
el Οδηγίες χρήσης - Συμπιεστικό ένδυμα για τη θεραπεία οιδημάτων
cs Návod k použití - Kompresivní prádlo pro terapii edému
pl Sposób użytkowania - Odzież kompresyjna do leczenia obrzęków
et Kasutusjuhend - Kompressioonriietus turse raviks
lv Lietošanas instrukcija - Kompresijas apģērbs tūskas ārstēšanai
lt Naudojimo instrukcija - Kompresiniai drabužiai edemai gydyti
sk Návod na použitie - Kompresívne odevy na liečbu opuchov
sl Navodila za uporabo - Kompresijsko oblačilo za terapijo za edeme
hu Használati utasítás - Ödématerápiás kompressziós ruházat
ro Instrucțiuni de utilizare - Îmbrăcăminte compresivă pentru edemoterapie
bg Инструкция за употреба - Компresионно облекло за терапия на отоци
hr Upute za uporabu - Kompresijska odjeća za terapiju edema
ru Инструкция по использованию - Изделия из компрессионного трикотажа для лечения отеков
tr Kullanım talimatı - Ödem tedavisi için kompresyon kıyafeti
ko 독일어 - 부종 치료를 위한 압박복
ja 使用説明書 - 水腫治療用コンプレッションガーメント
ar مروت لالجالعال طغاض سبالم - لامعت سالال ااداشرا

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Compression garments for oedema therapy

en

Purpose

Compression stockings are the basic therapy in the treatment of venous disorders and oedemas. In addition to their medical effect, compression garments support your well-being as a whole. Compression therapy comprises measures of exterior, planar pressure applications in case of venous disorders and lymphatic drainage problems. The application of planar pressure is to prevent the development of oedemas and support venous reflux or lymphatic drainage. Compression therapy generally is a long-term treatment and requires the corresponding products which – if the medical indication is given – the patient him/herself applies upon the doctor's instruction.

Indications

Primary and secondary varicosis; varices during pregnancy; supporting sclerotherapy; after surgical venous interventions; superficial thrombophlebitis, as well as condition after healed up phlebitis; deep vein thrombosis; condition after thrombosis; postthrombotic syndrome; thrombosis prophylaxis for mobile patients; CVI of stages I – III according to Widmer and C1 – C6 according to CEAP; ulcer prevention; venous insufficiency; angiodysplasias; lymphoedemas; oedemas during pregnancy; post-traumatic oedemas; postoperative oedemas; cyclic idiopathic oedemas; lipoedemas; lipo-lymphoedemas; phlebo-lymphatic oedemas; congestive conditions as a result of immobility (arthrogenous congestion syndrome; pareses and partial pareses of extremities); improvement of venous symptoms; occupational oedemas (walking and standing work; sedentary work).

After burns and for scar treatment, we recommend our special products for scar therapy.

Required combination supplies:

- Compression Capri pants with compression calf stocking
- Compression Bermuda pants with compression half-thigh stocking
- Compression leggings with compression calf stocking
- Compression cycling pants with compression thigh stocking

Absolute contraindications

Peripheral arterial occlusive disease (PAOD) stage III/IV; decompensated heart failure; septic phlebitis; phlegmasia coerulea dolens; acute erysipelas; acute bacterial, viral or allergic inflammations; untreated or advanced dermatoses /contact dermatitis; untreated open wounds; swelling of extremities of unknown origin; untreated malignant tumours.

Do not use Silon-TEX® in case of fissured skin or open wounds; also, do not use it for transplants which did not heal in a stable manner.

Relative contraindications

Pronounced exudative dermatoses; malignant lymph oedema; peripheral arterial occlusive disease (PAOD) stage I/II; peripheral neuropathy associated with diabetes mellitus; sensitivity disorders of extremities; compensated heart failure; skin infections; dermatoses (e.g. blistering dermatoses); intolerance of pressure or of product ingredients; primary chronic polyarthritis. In case of open wounds in the treatment area, the wound area must first be covered with an adequate wound dressing or a bandage before the compression garment is applied.

Thuasne Deutschland GmbH accepts no liability for failure to observe these contraindications.

Side effects

No side effects are known if properly used. In case of any negative changes during the application, please consult your doctor or the medical retailer immediately. In case of complaints in connection with the product, such as damage of the knitted fabric or lack of fit, please contact directly your medical retailer. Please note the following if Silon-TEX® is sewn into your compression garment: Should negative changes occur, such as itching or a skin allergy due to silicone intolerance, please contact immediately your doctor or your medical retailer. Should incompatibilities regarding one or several ingredients of this product be known, please consult your doctor prior to use. If skin irritations or allergies persist, stop using Silon-TEX®.

The manufacturer shall not be liable for damages / injuries which resulted due to improper handling or misuse.

Precautions

Follow the instructions given by your doctor, by an orthopaedic technician or medical expert personnel. If you experience severe pain or discomfort, please see your doctor, an orthopaedic technician or medical expert personnel. Store at room temperature, preferably in the original packaging.

To guarantee its effectiveness and for reasons of hygiene, the product may not be passed on to other patients. Any serious incident that occurs in connection with the product must be reported to the man-

ufacturer and to the competent authority of the Member State in which the user and/or patient are/is domiciled.

Directions for putting the garment on

Your medical retailer will be pleased to show you the correct application.

Rough, unfiled fingernails or sharp-edged rings etc. may damage the extremely fine knitted fabric. Therefore wear rubber gloves please. We recommend our fitting aids if there are difficulties putting on the garment.

If Silon-TEX® is sewn into your compression garment, please note the following:

- Scar and the surrounding skin must be thoroughly cleaned.
- The scar tissue is to be kept dry and grease-free.
- Do not use ointments, creams, gels or lotions in combination with Silon-TEX®.
- Put on your compression garment such that the Silon-TEX® is placed directly on the scar area and fits smoothly.
- Do not use in case of skin irritations.
- Do not use in case of a known intolerance of the ingredient (high-purity silicone).
- Should skin incompatibilities occur during application, please stop using the product and consult your doctor
- Do not use Silon-TEX® on open wounds!

Care

Compression garments should be washed daily, if possible. The garment fabric not only collects remainders of ointments and grease but also skin scales, sweat and dust. Regular washing will protect you from infections and your garment from wear. You need a spare set of garments so that you need not interrupt your compression therapy while the other set is washed and dried.

Always turn your compression garment inside out to wash it. Use a mild detergent for washing by hand or for machine wash (gentle cycle at 40° C) without optical brightener and fabric softener. In a laundry bag, compression garments are specially protected. In order to reduce the drying period, put the compression garments into a terry cloth towel, roll it up firmly and press the garment out (do not wring it out). Do not leave the moist compression garment in the wet towel. Hang up your garment loosely and leave it to air-dry. Never put it in the sun or onto a radiator to dry.

If Silon-TEX® is sewn into your compression garment, please dry the piece of garment especially well for reasons of hygiene and store it loosely as far as possible and with ventilation.

Do not expose compression garments to any damaging effects, such as: strong sunlight, intense heat, ointments, oils and greases, stain remover, gasoline and aggressive detergents.

Material composition

Please find the exact information on the sewn-in fabric label.

Note: The sewn-in fabric label is an important component of the compression garment. If you remove this labelling from a Thuasne product, any claim for warranty, repair or replacement shall be waived.

Material composition of Silon-TEX®: Silicone layer: 100 % medical grade silicone, carrier material: 83 % poly-amide, 17 % elasthan

Storage information and duration of use

Store in a dry place and protect from sunlight. For series products, a maximum shelf life of 36 months applies. This is followed by the duration of use of 6 months for the compression garment. The hourglass pictogram shows the duration of use printed on the carton label. Customised manufacture concerns products for immediate use.

For serial sizes and custom-tailored garments, we guarantee their medical efficacy for a maximum period of wear of six months. Prerequisite for this will be the proper treatment of the compression materials (for example, their care or putting them on or taking them off). If necessary, your doctor should then prescribe new compression garments. In case of any significant change of relevant body sizes (e.g. due to a successful therapy, weight change), new measurements and a new supply may also be justified earlier already.

Disposal

You may dispose of your medical compression garment via the residual or household waste. No special criteria exist for their elimination.

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- IL** THUASNE CERVITEX Ltd. · The Bell House, 20 Hataas St. · Kfar Saba, 44425
- S** THUASNE SCANDINAVIA AB · Box 95 · 131 07 Nacka
Thuasne UK Ltd
THUASNE DEUTSCHLAND GmbH