

THUASNE

XLR8

Botte de marche pneumatique à coque rigide, version courte ou version longue
Inflatable rigid-shell walking boot, short version or long version

- de** Unterschenkel-Fußorthese in Schalenform, mit hohem oder niedrigen Schaft und Luftpolstern
- nl** Pneumatische walker met rigide schaal. Verkrijgbaar in korte versie of lange versie
- it** Stivale da deambulazione pneumatico rigido, versione breve o versione lunga
- es** Bota de caminar neumática de carcasa rígida, versión corta o versión larga
- pt** Bota de caminhada pneumática de invólucro rígido, versão curta ou versão longa
- no** Airwalker (gå-ortose) med stiv ytterstøvel, kort eller lang utgave
- da** Luft-støttestøvle med stiv skal i kort eller lang udgave
- fi** Jäykkärunkoinen paineilmatäytteinen kävelysaapas, lyhyt tai pitkä malli
- sv** Pneumatiskt Walker fotstöd med styvt skal, kort eller lång version
- el** Νάρθηκας ποδοκνημικής πνευματικός, με άκαμπτο περίβλημα, κοντή έκδοση ή ψηλή έκδοση
- cs** Pneumatická fixační ortéza s pevnou výztuhou, krátká nebo dlouhá verze
- pl** Pneumatyczny but ortopedyczny do chodzenia ze sztywną skorupą, wersja krótka lub długa
- lv** Stingra apvalka piepūšams pastaigu zābaks – īsā vai garā versija
- lt** Pripučiamas eiti padedantis įtvaras standžiu korpusu, ilgas arba trumpas
- et** Jäiga korpusega suruõhksaabas, madal või kõrg variant
- sl** Škorenj opornica z napihljivo nogavico in togim ogrodjem, nizka ali visoka
- sk** Pneumatický pevný walker, krátka alebo dlhá verzia
- hu** Merev héjjal rendelkező pneumatikus járócsizma, rövid változat vagy hosszú változat
- bg** Пневматичен ботуш за ходене с твърд корпус, къс или дълъг вариант
- ro** Cizmă pneumatică cu carcasă rigidă, versiune scurtă sau lungă
- uk** Пневматичний ортез гомілковостопного суглобу з твердим каркасом, висока та низька модель
- zh** 充气式硬壳步行靴，短款和长款
- hr** Pneumatska "čizma" za hodanje s čvrstim kalupom, kratka ili duga
- ms** But berjalan kaki cangkerang tegar boleh kembang, versi pendek atau versi panjang
- ar** جزمة مَشِي هوائية ذات هيكل صلب، قصيرة أو طويلة

1

2

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Inflatable rigid-shell walking boot

Description:

- XLR8 is an inflatable rigid-shell walking boot. It consists of a main rigid plastic part, fastened using straps, and a protective foam bootee for the lower leg and foot.
- Inside the rigid shell, on the back part, there are two independent, reinforced air cells that can be inflated and deflated using two pumps (also independent) integrated on the outside of the boot, on the right and left.
- The rigid strap system integrated on the rigid part is used to tighten the device and hold it in place.
- Product can be worn on either foot.

Composition:

Rigid shell: Polypropylene
Bootee: Polyurethane
Straps: Polyamide

Mechanism of action/Purpose:

The walking boot is designed to be used on the lower limb by adults, in the indications outlined below. It is used to immobilize and stabilize the lower third of the leg, ankle and foot.

Properties:

- Full rigid shell walking boot to support and protect the lower third of the leg, ankle and foot.
- Two integrated independent air cells ensure an anatomical fit and aid oedema reduction.
- Curved, anti-skid sole and adjustable toe guard (3 possible positions).

Indications:

- Severe ankle sprains
- Stable fractures or internally fixed fractures of the lower third of the leg, ankle or foot.

The choice of the short or long version of the walking boot is at the discretion of the health professional.

Contraindications:

- Not suitable for unstable fractures or fractures of the upper part of the tibia or fibula.
- Do not use the product in the event of an uncertain diagnosis.
- Do not apply the product in direct contact with broken skin.
- Not suitable for patients weighing > 160 kg.
- Known allergy to any of the components.

Putting on the walking boot:

- Check the product is not damaged before use.
- Make sure the product is the right size to fit the patient.
- The systematic use of a long sock is recommended when wearing the walking boot. The sock should be worn against the skin, under the walking boot, and the top of the sock must be above the top of the boot, if possible.
- When using for the first time, make sure you remove the plastic elements placed over the inflator bulbs to protect them during transport.
- The health professional must supervise fitting of the product and the patient's specific walking conditions when the device is used for the first time.
- Before putting on the walking boot, make sure the tightening system has been fully unfastened and that the air cells are fully deflated.

1. Remove the walking boot from its plastic cover, unfasten the touch fastening straps, remove the rigid front part and open the foam bootee.
2. Sit down and carefully place your lower leg and foot in the foam bootee, checking that your heel is correctly positioned at the bottom of the boot (A). Close the foam bootee, applying it firmly around the lower leg and foot (B).
3. Put the rigid front part back in place (C).
4. Fasten the walking boot again, tightening the touch fastening straps from bottom (D) to top (E). To find the perfect level of tightness, insert the top touch fastening strap through the loop cut directly in the rigid shell (E1).
5. Adjust the toe guard.
6. To increase compression while maintaining comfort at the ankle bones, two air cells are integrated on either side of the boot, to be inflated by pressing the blue inflator bulb (F-F1). To deflate the cells (and reduce compression), open the cap on the back of the inflator bulb and press the valve to release air (G). Do not hesitate to press on the air cell to help release air.
7. Try walking with the boot to check that it fits properly and that it is working correctly.

Taking off the walking boot:

- Loosen the straps and remove the rigid front part.
- Open the bootee up completely to remove the lower leg, leaving the bootee inside the boot.

Precautions:

- Follow the advice of the health professional who prescribed or supplied the product.
- Do not apply the device in direct contact with a greasy substance (ointment, cream, etc.)
- It is recommended to adequately tighten the product to achieve immobilization without restricting blood. In the event of unusual sensations, loosen the straps and/or deflate the air cells and ask a health professional for advice.
- Check the condition of the leg and foot daily (paying particular attention to patients with loss of feeling in the leg). In the event of discomfort, numbness affecting the leg, unease, pain, swelling, changes in sensation, unusual reactions (allergies, etc.) or changes in skin colour, discontinue use of the device and consult a health professional.
- Driving while wearing XLR8 is not recommended. The patient's ability to drive while wearing the walking boot should be assessed by the health professional.
- If the device comes into contact with water, dry the fabric part and carefully wipe the rigid plastic part using a dry cloth.
- For hygiene and performance reasons, do not re-use the product for another patient.
- Any serious incidents related to the device should be reported to the manufacturer and to the competent authority of the member state in which the user and/or patient is resident.

Unwanted side effects:

- Possible risk of itching/skin reactions following contact by the product with skin
- Possible risk of venous thrombosis
- Possible risk of back pain due to gait pattern modification.

Maintenance / Washing instructions:

- Product can be washed by following the instructions in this leaflet and on the labelling.
- Store at room temperature, preferably in its original packaging.
- Foam bootee fully removable for washing.
- Hand wash the bootee in cold water.
- Do not use detergents, fabric softeners or abrasive products (products containing chlorine).
- Squeeze out excess water. Air dry away from heat sources.
- Wash the rigid shells with a damp sponge.

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Thuasne Deutschland GmbH - Im Steinkamp 12 - 30938 Burgwedel - Deutschland
Thuasne Benelux B.V. - Klompenmaker 1-3, 3861 SK Nijkerk - Nederland
Thuasne Scandinavia AB - Ryssviksvägen 2, 131 36 Nacka - Sweden
Thuasne Cervitex Ltd - 6 Hanagar St., POB 2447, Kfar Saba 4442536 - Israel
Thuasne UK Ltd - Unit 4 Orchard Business Centre - North Farm Road - Tunbridge Wells - Kent TN2 3XF - United Kingdom