



RETAILER IDENTIFICATION

Customer code:

Case No. for renewal

Patient's surname:


Patient's first name : 1st treatment

Gender : M F Child

Date:

Patient's height : Quantity:

If possible, please enclose photos of the limb to be fitted.
Please draw in the contours of the garment on the diagram and cross unnecessary measures.
Please also fill in torso form No. 1b.



circumference in cm
 length in cm
 half-circumference in cm

35 and 36 : For the body only
 Measurement taken at an angle for a brief-cut body
 Measurement taken straight for a short-cut body

Please contact your regular Thuasne distributor



RETAILER IDENTIFICATION

Case No. for renewal

Patient's surname:

Patient's first name : 1st treatment

Gender : M F Child Date:

Customer code: Patient's height : Quantity:

If possible, please enclose photos of the limb to be fitted.
Please draw in the contours of the garment on the diagram and cross unnecessary measures.
Please also fill in torso form No. 1a.

Indications

- Major burn
 EDS (Ehlers-Danlos Syndrome)

Models

- Shoulder pad
 Body
 Vest

Fabrics

- Cicatrex Nature
 Cicatrex Filifine
 Cicatrex AirSkin

	WHITE	BEIGE	BLACK
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>

Options

Center opening

- Zip
 Front position Back position
 Hook to help fasten the vest (by default)

Armpit compression

- Right Yes No*
 Left Yes No*
 * Soft fabric by default

Collar

- Round neck
 "V" neck
 Front collar cut height (from the sternum): cm
 Back collar cut height (from the 7th cervical vertebra): cm
 Stand-up collar
 Circumference at top: cm
 Circumference at bottom: cm
 Front height: cm
 Back height: cm
 Top of sternum/neck height: cm

Cups

- With compression
 Without compression
 Depth: A B C D E F

Velcro fastening (for panty)

- Vest over panty
 Panty over vest

Additional opening**

- Position: (to be drawn in)
 Zip Velcro Length: cm

Soft fabric**

- Elbow
 Other

Other**

- Compression pad
 Position: (to be drawn in)
 Dimensions: cm
 Belt

** Extras, not included in the cost

RIGHT SIDE

13 Shoulder length from the acromion process to neck circumference No. 1

14 Armhole circumference

16 Top of sleeve

15 Arm circumference at the armpit

18 Length to the middle of the upper arm

17 Upper arm circumference (widest part)

20 Length to the elbow

19 Elbow circumference

22 Length to the middle of the forearm

21 Forearm circumference (widest part)

24 Length to the wrist

23 Wrist circumference

36 Thigh circumference*

1 Neck circumference from the 7th cervical vertebra to the top of the sternum

2 Chest width, breathing in (armpit to armpit)

4 Top of sternum - armpit line length

3 Chest perimeter, breathing in

6 Top of sternum - under the bust length

5 Circumference under the bust, breathing in

8 Top of sternum - waist length

7 Waist circumference

10 Top of sternum - bottom of vest length (for a body, length to the crotch)

9 Hip circumference

37 Back length (7th cervical vertebra - bottom of vest). For a body, length to the crotch, measure with the tape measure flat against the skin

39 Circumference tip of the breast breathing in

38 Length top of sternum - tip of the breast

a** Length neck base - tip of the breast

b** Length neck base - under the bust

LEFT SIDE

12 Armhole circumference

11 Shoulder length from the acromion process to neck circumference No. 1

25 Arm circumference at the armpit

26 Top of sleeve

27 Upper arm circumference (widest part)

28 Length to the middle of the upper arm

29 Elbow circumference

30 Length to the elbow

31 Forearm circumference (widest part)

32 Length to the middle of the forearm

33 Wrist circumference

34 Length to the wrist

35 Thigh circumference*

* For the body only
 ** To fill out only from D cup

Please contact your regular Thuasne distributor