



RETAILER IDENTIFICATION

Customer code:

Case No. for renewal

Patient's surname:

Patient's first name : 1st treatment

Gender : M F Child Date:

Patient's height : Quantity:

**If possible, please enclose photos of the limb to be fitted.
 Please draw in the contours of the garment on the diagram and cross unnecessary measures.
 For open fingers, please indicate the desired length for each finger.**

RIGHT HAND LEFT HAND
 Fill in one form for each side

Indications

- Major burn
- EDS (Ehlers-Danlos Syndrome)

Fabrics

- | | | | |
|---|--------------------------|--------------------------|--------------------------|
| | WHITE | BEIGE | BLACK |
| <input type="checkbox"/> Cicatrex Filifine* | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| <input type="checkbox"/> Cicatrex AirSkin** | - | <input type="checkbox"/> | <input type="checkbox"/> |

*Minimum circumferences of 4.6 cm, minimum length of fingers of 2 cm.
 **Minimum circumferences of 2.6 cm, minimum length of fingers of 1 cm.

Options

Proximal end

- With length on wrist (5 cm max)
- With sleeve beyond 5 cm and up to the armpit as maximum (fill in form ARM No. 2) Separate Attached

Finger protector

- Finger protector with strap
- Finger protector without strap

Opening


- Type: Zip Velcro Dimensions: cm
- Position: Palm of hand Back of hand Little finger side

Other

- Compression pad
- Dimension (l x h): cm
- Position: to be drawn in Removable Fixed
- Palm of hand Back of hand
- Palm support

Comments

Circumferences		Measures in cm
22	wrist circumference + 5 cm max	
23	wrist circumference	
E	circumference of the hand at the head of the metacarpal bones	
F	circumference of the little finger at the 1st phalanx	
G	circumference of the little finger at the 3rd phalanx	
H	circumference of the ring finger at the 1st phalanx	
I	circumference of the ring finger at the 3rd phalanx	
J	circumference of the middle finger at the 1st phalanx	
K	circumference of the middle finger at the 3rd phalanx	
L	circumference of the index finger at the 1st phalanx	
M	circumference of the index finger at the 3rd phalanx	
N	circumference of the thumb at the 1st phalanx	
O	circumference of the thumb at the 2nd phalanx	
Desired length	Fingertips	Measures in cm
D	Length of the palm surface of the hand	
P	Length of the little finger	<input type="checkbox"/> OPEN <input type="checkbox"/> CLOSED
Q	Length of the ring finger	<input type="checkbox"/> OPEN <input type="checkbox"/> CLOSED
R	Length of the middle finger	<input type="checkbox"/> OPEN <input type="checkbox"/> CLOSED
S	Length of the index finger	<input type="checkbox"/> OPEN <input type="checkbox"/> CLOSED
T	Length from the tip of the index finger to the join with the thumb	
U	Thumb length	<input type="checkbox"/> OPEN <input type="checkbox"/> CLOSED
Z	Length of the sleeve up to 5 cm max	



circumference in cm length in cm

Please contact your regular Thuasne distributor